



INTELLIGENT CONTROLLER DUAL SUSPENSION SYSTEM

WATTI-YE-246



BRUSHLESS LIGHTWEIGHT POWER WHEELCHAIR

This Power wheelchair is a most light-weight one with a smallest folding size. the chair weights 21.6 kg and with a max loading capacity of 135 kg. It has the best Brushless motor and Intelligent controller. The enforced frame and bigger seat and higher backrest improve the comfort and increase the max loading capacity. A patented folding system makes the FOLD & UNFOLD in 3 steps, in 1 second.

FEATURES

- Top quality Brushless motors, quiet and drive longer
- Awesome Suspension system offers super comfort
- 21.6kg only, easy and effortless transport and fit in all vehicles
- Brushless controller, intelligent control and programming
- Li-ion battery from Japn, max driving range to be 18km
- Easy fold and go; flight free for travel

YATTLI INDUSTRY CO., LTD.

Tel: +86 755 8170 7199

Web: www.yattli.com

Email: info@yattli.com

Factory: No.2 Lingchang Road, PingShan Village, Tang Xia Town, DongGuang City, GuangDong, China

Office: RM.1403, BLK.B of Duocai Office Complex, Guanlan High-tech Zone, Longhua DIST., Shenzhen, China



USED IN VARIOUS
LIFE SCENES



Product Details



Brushless controller



Adjustable armrest



Top quality brushless



Suspension system



Packing Bag



YATTLL INDUSTRY CO., LTD.

Tel: +86 755 8170 7199

Web: www.yattll.com

Email: info@yattll.com

Factory: No.2 Lingchang Road, PingShan Village, Tang Xia Town, DongGuang City, GuangDong, China

Office: RM.1403, BLK.B of Duocai Office Complex, Guanlan High-tech Zone, Longhua DIST., Shenzhen, China



Footnotes

- 1) Range and speed vary with userweight-terrain type, battery charge battery condition and tyre pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- 3) Battery weight may vary based on manufacturer.

Parameter

Front wheel: 7", PU Solid Tyre

Rear Wheel: 8", PU Solid Tyre

Net Weight: 21.6kg

Brake: Electromagnetic Braking system

Batteries: Lithium Battery 24V 11Ah

Motors (pair): Brushless Motor 250W

Charger: High Power 24V 2A

Controller: Reliable Brushless Controller

Maximum Range: 15-18 km

Suspension System: Yes

Ground Clearance: 8 cm

Max Speed: 8 km/h

Charging Hour: 5-6 Hours

Max Slope: 12 degrees

Weight Capacity: 135 KG

Dimension (Folded): 71*58*30.5cm

Dimension (Folding): 99*58*98cm